

Looking for a way to Manage, Train, Certify, and Record complex processes?

Assembly and Inspection Tracker

(See our website for details)

Focal*Point Release v. 4.14.8

Features and Enhancements

- **Expanded Wireless Compatibility** – Updated COMM engine for expanded device compatibility.
- **Enhancements Include** – Port re-connect maintains FTP and Timer folder status, updated RFA file overwrite processing, corrects errors in the file revision process, eliminates machine status display errors. Corrects the default socket / port numbers for Win7 compatibility. Includes updated reports.
- **RFA** – Updated to provide better file matching when selecting right and left pane program items.
- **FTP and Data Collection** – Allows serial DPRNT data capture when using FTP for program transfers.
- **Timer Folders** – Provides a managed approach to network shared folders. Includes program redirect for RS232 connected machines. Excellent for use with Mazak CMT remote file requests.
- **FTP** – Focal*Points FTP server supports Managed, Restricted and Unrestricted operations.
- **Numeric DNC ID Preference** – The Settings menu now has an option to specify the “highest” numeric value or to “re-use” values that have been deleted.
- **DNC Library Functions** – The Add Entry form (RFA Included) has been completely redesigned to include an embedded program viewer, DNC ID comments with enhanced ID options.
- **New Windows 7 Compatibility** – Includes expanded help for setup and operation with Win7.
- **Multi-User Administration has been completely revised** - New options provide better control in a multi-user environment including a method to re-assign a Cell ID when a PC is rebuilt or replaced.
- **New GForce Options** - The GForce editor will open several programs from the DNC Library, and provides options to print information headers.
- **Received File Administration Enhancements** - The RFA will now display orphaned files and retains last user information.
- **New Data Collection (Option)** – Capture DPRNT statements containing macro and probe data. Append date and time stamps for calculating cycle times.
- **Enhanced Auto Library Transfer** – Added features to check for locked files and tag particular extensions for deletion or normal processing as well moving files to a relocation folder.

For more information call 508-523-1763 or e-mail inquiries to info@cncinnov.com.

Increase Productivity - Ask About our (OEE) Machine Asset Monitor